



Automation and Process Control for Ultra High Throughput Screening

The Automation Partnership designs and manufactures advanced automation systems for research, development and manufacturing processes within the life science industries.

Business Problem

High Throughput Screening (HTS) is an important part of the drug discovery process. This allows a company to screen the initial effect of thousands of components against identified targets ('diseases'). Many pharmaceutical companies have set minimum goals of running 100s of targets per year against 500,000 compounds per target. To reach this goal, an average of 200,000 compounds must be screened every working day of the year.

Until recently, the only choice available to HTS management to meet these goals was to buy more laboratory robotics, hire more skilled staff to support the robotics and find the additional office and lab space to house both. This option is both difficult and costly.

Tessella Solution

The Automation Partnership, in close collaboration with AstraZeneca (Alderley Park, UK), have developed Asset, a fully integrated, HTS system that can process a plate every 30 seconds, using any assay technology. Tessella, with their background of involvement in the pharmaceutical industry, has helped the Automation Partnership to develop the software required to control Asset.

Complex assays that require multiple plate washes, incubations, reagent additions and even multiple measurements are easily set-up in the user defined assay protocols.

The control software then takes over, coordinating the movement of each of the robots within the system, monitoring the environmental sensors, and activating the plate readers. In addition, maintenance panels for each of the robots allow engineers to test the robots and adjust them during and after maintenance.



The software was developed using Visual C++. It runs on a bank of computers running Windows NT, which communicate with each other using DCOM (Distributed Component Object Model).

Results and Benefits

Asset enables more screens to be completed faster with a greater number of compounds than is possible with previous approaches. Reduced operator involvement and closely controlled processing eliminates the variability problems associated with batch processing. The reliability of both

the hardware and the software mean that Asset can be left unattended for overnight operation, further increasing the throughput that can be achieved.

Tessella plc 26 The Quadrant, Abingdon Science Park, Abingdon, Oxfordshire OX14 3YS, UK
T: +44 (0)1235 555511 | F: +44 (0)1235 553301 | E: info@tessella.com

Tessella Inc 233 Needham Street, Suite 300, Newton, MA 02464, USA
T: 1 617 454 1220 | F: 1 617 454 1001 | E: info@tessella.com

Tessella – successfully delivering IT and consulting services to world leaders in R&D, science and engineering.

For decades, Tessella has been successfully delivering IT and consulting services to world leaders in R&D, science, and engineering. Through the application of scientific methods and rigorous quality procedures, we enable clients in life sciences, energy, the public sector, and consumer industries to achieve a wide range of objectives, including, forecasting floods, developing fusion power, enhancing military sensor capability, increasing drug discovery and development efficiency, and reducing risk to health and the environment in the extraction and production of oil and gas. With offices in Europe and North America, global companies rely on Tessella for business critical assignments.

Copyright © Tessella plc 2009, all trademarks acknowledged.

